

SPECIAL IMPLANT RESTORATION RX OFFER**

1. Send this Rx with your next crown over implant case to save \$30 per crown or implant abutment (up to 3 units).
2. Carefully package your case, including this Rx, and tape box securely closed.
3. To schedule shipping pickup, call us at 888-336-1301.

Offer expires March 31, 2022. **Price does not include round-trip shipping or applicable taxes. Limit one specially priced Rx per case. Web offer can only be used a maximum of four times per account. Special pricing not valid with any other offer.

Dr. Name _____ Acct. # _____

Patient ID/Name _____ Male Female
First Last

Deliver by 5 p.m. on _____ Call before starting case


WEB RX

\$30 off Crown Over Implant

\$30 off Implant Abutment



FINAL SHADE INSTRUCTIONS



Tooth No. _____

Final Shade _____

Occlusal Staining: Light* Med Dark None

SELECT RESTORATION TYPE

Crown Splinted Crowns Bridge

ALL-CERAMIC RESTORATIONS

- BruxZir Full-Strength* (1,150 MPa) **NEW!** BruxZir Esthetic (870 MPa) (stump shade recommended for restorations less than 1.5 mm thick)
- PrismaTik Clinical Zirconia IPS e.max

SELECT ABUTMENT TYPE

Custom Implant Abutments

- Titanium Abutment* Zirconia w/ Ti-Base
- Prepare existing abutment
- See reverse for implant systems supported


SCREW-RETAINED RESTORATIONS

- BruxZir Full-Strength* (1,150 MPa) (w/ Ti-Base) IPS e.max
- NEW!** BruxZir Esthetic (870 MPa) (w/ Ti-Base)

PARALLEL ABUTMENTS

- No Yes (Indicate which abutments will have restorations splinted together for insertion.)

ABUTMENT MARGIN DEPTH



If left blank, default values will be used.


CONTOUR AND OCCLUSION DESIGN

- Embrasures: Open Closed*
- Occlusion: Heavy Light*
- Contacts: Pinpoint Light Broad & Tight*


RESTORATION PONTIC DESIGN



RESTORATION MARGIN DESIGN



METAL DESIGN



IF NO OCCLUSAL CLEARANCE

- Call doctor Spot opposing
- Make this a permanent note in my master file

*Standard unless specified otherwise


Signature _____

License # _____ Date _____


Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details.

Indicate implant system _____


Indicate diameter below



ABUTMENT MARGIN DESIGN



ABUTMENT EMERGENCE PROFILE



TERMS AND WARRANTY INFORMATION

We honor VISA, MASTERCARD, AMEX and DISCOVER.

TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance. Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: DigiTech Dental Restorations is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

LIMITED WARRANTY/LIMITATION OF LIABILITY. For warranty terms and conditions and limitation of liability, visit digitechlab.com/privacy-policy-warranty.



• BruxZir® Restorations




• Custom Abutments



• All-Ceramic Restorations
• PFM Restorations

PREPARATION GUIDELINES



BruxZir Esthetic

- A. 1.25 mm ideal reduction (0.7 mm minimum)
- B. Chamfer or modified shoulder margins preferred
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

BruxZir Full-Strength

- A. 1.0 mm ideal reduction (0.5 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

FLAT-RATE PRICE ON CUSTOM ABUTMENTS AND SCREW-RETAINED CROWNS IS AVAILABLE FOR THE FOLLOWING IMPLANT SYSTEMS

BIOMET 3i™
Certain®

CAMLOG®
SCREW-LINE

DENTSPLY Implants
ANKYLOS® C/X
ASTRA TECH Implant System®
ASTRA TECH Implant System® EV

Glidewell Direct
Hahn™ Tapered Implant System
Inclusive™ Tapered Implant System

HIOSEN®
HG System

MegaGen
AnyRidge® Implant System

Nobel Biocare
Brånemark System® RP
NobelActive®
NobelReplace®

Straumann®
Bone Level
Tissue Level

Zimmer Dental
Screw-Vent®

Restorations for the implant systems above are fabricated using components manufactured by Prismatic Dentalcraft, Inc. Prices may vary for other implant systems due to the need for original equipment manufacturer (OEM) components. Inclusive is a registered trademark of Prismatic Dentalcraft, Inc. Hahn Tapered Implant is a trademark of Prismatic Dentalcraft, Inc. All other trademarks are property of their respective owners.

All rush cases must be prescheduled by calling **888-336-1301** before the case is shipped. Time of pickup and delivery may affect turnaround time.

