

Dr. Name \_\_\_\_\_ Phone # \_\_\_\_\_

Account # \_\_\_\_\_ Patient ID/Name \_\_\_\_\_

First

Last

Address \_\_\_\_\_

Email \_\_\_\_\_ Deliver by 5 p.m. on \_\_\_\_\_

Enclosed with Case:  Impressions  Models  Bite  Photos  Other: \_\_\_\_\_




Indicate implant system \_\_\_\_\_  
See reverse for flat-rate pricing info

Indicate implant diameter \_\_\_\_\_ mm

Please order all implant components for this case



### SHADE INSTRUCTIONS



Tooth No. \_\_\_\_\_

Final Shade \_\_\_\_\_

Occlusal Staining:  Light\*  Med  Dark  None

### SELECT FINAL CUSTOM ABUTMENT

Choose Material

- Titanium\*
- Gold-Colored Titanium
- Zirconia w/Ti-Base
- Gold Alloy

Choose Abutment


- Glidewell\*
- BIOMET 3i Encode

Prepare existing abutment

### PARALLEL ABUTMENTS


No  Yes (indicate which abutments will have restorations splinted together for insertion)

### ABUTMENT MARGIN DEPTH



If left blank, default values will be used

### ABUTMENT EMERGENCE PROFILE



Surgical placement  Tissue displacement\*  No tissue displacement

### BIOTEMPS PROVISIONALS

- Cementable
- Screw-Retained
- Reinforcement:  None  Wire  Fiber  Cast-Metal

### SELECT RESTORATION TYPE

Crown  Splinted Crowns  Bridge

### CEMENT-RETAINED RESTORATIONS

- BruxZir Full-Strength\* (> 1,000 MPa)
- BruxZir Esthetic (870 MPa)
- Bilayered Clinical Zirconia
- IPS e.max (anterior)
- IPS e.max (posterior)
- Obsidian All-Ceramic
- Obsidian Fused to Non-Precious
- Obsidian Fused to Noble
- Obsidian Fused to White High Noble

### SCREW-RETAINED RESTORATIONS

- BruxZir Full-Strength\* (w/Ti-Base)
- BruxZir Esthetic (w/Ti-Base)
- IPS e.max (w/Ti-Base)
- Obsidian All-Ceramic (w/Ti-Base)
- Obsidian Fused to Noble
- Obsidian Fused to White High Noble
- Bilayered Clinical Zirconia (w/Ti-Base)

### CONTOUR AND OCCLUSION DESIGN

Embrasures:  Closed\*  Open  
Occlusion:  Light\*  Ideal  Open \_\_\_\_ mm  Out  
Contacts:  Broad & Tight\*  Pinpoint  Light

### PONTIC DESIGN



### ESSIX RETAINER

Scan & Save (Fee Applies)  
Thickness:  0.8 mm  1.0 mm\*  1.5 mm  
Essix Design:  Standard  Scalloped  
Palatal Coverage:  Yes  No

### IF NO OCCLUSAL CLEARANCE

Call doctor  Spot opposing  Metal occlusion  
 Metal island  Make this a permanent note in my master file

\*Standard unless specified otherwise

Signature \_\_\_\_\_

Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details.

License # \_\_\_\_\_ Date \_\_\_\_\_

## TERMS AND WARRANTY INFORMATION

**We honor VISA, MASTERCARD, AMEX and DISCOVER.**

**TERMS:** Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance. Prices subject to change without notice. Rx must be enclosed with original case submission.

**NO-FAULT REMAKE POLICY:** DigiTech Dental Restorations is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

**LIMITED WARRANTY/LIMITATION OF LIABILITY.** For warranty terms and conditions and limitation of liability, visit [digitechlab.com/privacy-policy-warranty](http://digitechlab.com/privacy-policy-warranty).



• Custom Abutments



• BruxZir® Restorations



• All-Ceramic Restorations  
• PFM Restorations

### FLAT-RATE PRICE ON CUSTOM ABUTMENTS AND SCREW-RETAINED CROWNS IS AVAILABLE FOR THE FOLLOWING IMPLANT SYSTEMS

BioHorizons®  
Tapered Internal

BIOMET 3i™  
Certain®

CAMLOG®  
SCREW-LINE

DENTSPLY Implants  
ANKYLOS® C/X  
ASTRA TECH Implant System®  
ASTRA TECH Implant System® EV

Glidewell Direct  
Glidewell HT™ Implant System  
Hahn™ Tapered Implant System  
Inclusive® Tapered Implant System

HIOSSEN®  
HG System

MegaGen  
AnyRidge® Implant System

MIS®  
C1 Implant System

Nobel Biocare  
Brånemark System® RP  
NobelActive®  
NobelReplace®

Straumann®  
Bone Level  
Bone Level SC  
Tissue Level

Zimmer Dental  
Screw-Vent®

### FLAT-RATE PRICE ON CAD/CAM MILLED IMPLANT BARS IS AVAILABLE FOR THE FOLLOWING IMPLANT SYSTEMS

BIOMET 3i™  
Certain®

DENTSPLY Implants  
ASTRA TECH Implant System®

Nobel Biocare  
Brånemark System® RP  
NobelActive®  
NobelReplace®

Straumann®  
Bone Level

Zimmer Dental  
Screw-Vent®

Restorations for the implant systems above are fabricated using components manufactured by PrismaTik Dentalcraft, Inc. Prices may vary for other implant systems due to the need for original equipment manufacturer (OEM) components. Inclusive is a registered trademark of PrismaTik Dentalcraft, Inc. Glidewell HT and Hahn are trademarks of PrismaTik Dentalcraft, Inc. All other trademarks are property of their respective owners.

**All rush cases must be prescheduled** by calling **888-336-1301** before the case is shipped.  
Time of pickup and delivery may affect turnaround time.